



Measure for sizing



Two measurements to get the right size for your dog:

Measurement 1: Measure the circumference of your dog's leg above the knee just below the hip. This measurement should be at the highest point on the leg where you can measure all the way around the leg.

Measurement 2: Measure the circumference of the leg one inch above the hock (ankle) joint.

- **Medium:** Measurement 1: 11.25" Measurement 2: 4.25"
- **Large:** Measurement 1: 12.5" Measurement 2: 5.25"
- **Extra Large:** Measurement 1: 14" Measurement 2: 6"

*If your dog is in-between sizes, we recommend ordering the larger of the two sizes. It's okay if the brace overlaps slightly as long as the metal supports are on either side of the leg.

